



enhanced  
performance

Ref. No.  
**OPJ-00101-00**



TECHNICAL DATASHEET

## Lower Upper Set

Top and Bottom Sets used by Surgeons produced by Opia Medical are processed in the latest system and sterile environments. This product is an important part of PPE in healthcare services and is used to fulfill a dual function of preventing the contamination of microorganisms and body fluids from surgeons and operating room personnel to patients and vice versa. It is generally used before surgical operations. This allows Surgeons and Surgical nurses to move more comfortably during the operation.

### Technical Standards

It is breathable, waterproof, lightweight and easy to wear. Knitted or elastic cuffs are optionally used. Our products can be easily used by patient relatives, patients, nurses, doctors and laboratory workers. It consists of two parts, upper and lower. Optionally, the upper part consists of 2 or 3 pockets in one piece. The trousers part is desired with pockets and elastic. These features may vary according to the needs of our users.

Our products are produced as Sterile and Non Sterile.



Top Jersey

Size	Length	Width
S	68	50
M	72	58
L	75	65
XL	78	68
XXL	82	74

Lower Jersey

Size	Length	Width
S	102	48
M	106	53
L	108	55
XL	110	60
XXL	113	63

### TECHNICAL INFORMATION

ORDER REF. #	OPJ-00101-00
WEIGHT	40 Gsm
COLOURS AVAILABLE:	Medical Blue
SIZES AVAILABLE:	S, M, L, XL, XXL,
TYPICAL APPLICATIONS:	HOSPITALS, CONSTRUCTION INDUSTRY MANUFACTURE, FOOD PROCESSING, CHEMICAL MANUFACTURE
MATERIAL:	SMS Non-Woven
PER CASE:	80 PIECES (INDIVIDUAL POLYBAG)

EU-UDI COD Barcode



Product Barcode



Box Size



W: 400 L: 600 H: 400

Pallet: 80\*120 cm



Pallet: 24 Box

### TECHNICAL INFORMATION:



MDR 2696



### CERTIFICATION:

This product has been certified by: UDEM Adriatic d.o.o. Radnicka cesta 54/R3 10000 Zagreb, CROATIA  
Notified body: 2696 in accordance with Regulation (EU) 2017/745 of the European Parliament and of the council of 04 April 2022, and has been shown to comply with the standards: Centrocot ISO 13795-1  
ISO 22612-4, ATCC 9372, Ecoteks : AAMI PB-70, AATCC-42, AATCC-127, LEVEL II- III - IV

»» enhanced performance



Phone: +90 (850) 308 6742  
Email : sales@opiamedical.com  
Web : opiamedical.com



© 2021. Opia. Head Office: Opia Inovasyon Sağlık Medical Ar-Ge Hiz. Dan. San. Tic. A.Ş.  
| VAT Registration No. Gazikent 644 080 2138 | Registered in Gaziantep No. 62380  
İncilipinar Mah. 36006 Street. No: 20/29 | Şehitkamil | Gaziantep | TURKEY