



DISPOSABLE END TYPE 3/4 COVERALLS



Product Description

Opia Medical & Health, our disposable Industrial Type 3 and Type 4 overalls have a microporous structure that is extremely resistant to chemical particles and liquid chemical splashes. Our Laminated microporous Nonwoven fabrics and Industrial Type 3 and Type 4 overalls have a breathable structure. It is extremely comfortable for long-term use. Our processes, which start with the technical tests of the raw materials before production, are carried out with the label of quality control on the packaging by making the final quality control. All necessary standards are applied through the Occupational and Occupational Health, WHO (World Health Organization) and PPE Regulations in the processes up to the final stage of our industrial TYPE 3 and TYPE 4 overalls. The most important part of our products, which go through many stages in production, is keeping the sealing at the top level. Our products, especially used in Aspest applications, have over 98% protection. The combination of microporous Coverall, soft polypropylene non-woven fabric and microporous film allows moisture vapor to escape to keep the wearer comfortable. It has good barrier property for wet or liquid and dry particles. Ideal for Security, Mining, Clean room, Food industry, Medical, Laboratory, Pharmaceutical, Industrial pest control, Machine Maintenance and Agriculture. In addition, each of our products that are produced are checked one by one by the quality control unit and packaged with the "quality controlled" stamp.

Our products have MOR Certificate and GMP Certificates. In addition, it is produced according to the national medical production license.

Technical Standards

Material : SMS + Polypropylene+ Microporous film coated

Appearance: Hooded, elasticated wrists, waist and ankles.

Strength : Excellent resistance to liquid and chemical splash Impermeability:

All parts that are stitched are closed with a strip by applying a special application at certain temperatures.

Front Closure: Zipper system and on it, zipper closure according to the size of each product

SIZE	LENGTH	WIDTH	ARM LENGTH
S	185	155	58
M	190	160	60
L	193	163	63
XL	200	167	66
2XL	203	170	68
3XL	207	173	70



MDR CE 2696

CERTIFICATION NO: M.2022.MDR.1002

THIS PRODUCT:
 UDEM Adriatic D.O.O RADNICKA CESTA 54/R3 10000 ZAGREB CROATIA
 NOTIFICATION BOARD: 2696 CERTIFIED IN APRIL 2022 ACCORDING TO
 THE 2017/745 REGULATION OF THE COUNCIL OF THE EUROPEAN
 PARLIAMENT. COMPLIED WITH THESE STANDARDS:
 CENTROCOT ISO 13795-1, ISO 22612-4, AATCC 9372,
 EKOTEX : AAMI PB-70 ,AATCC -42 ,AATCC-127,LEVEL II-III-IV

ORDER REF. #	END TYPE 3-601 / END TYPE 4-602
WEIGHT	55-60 GSM
COLOURS AVAILABLE:	White, Red, Blue, Yellow
SIZES AVAILABLE:	S, M, L, XL, 2XL, 3XL
TYPICAL APPLICATIONS:	In Hospitals, Factories, Chemical Manufacturing
MATERIAL:	SMMS + Polypropylene + Microporous film coated
PER CASE:	50 PIECES
BOX SIZE:	40X60X40
PALLET SIZE:	80X120
PALLET BOX:	24



PHONE : +90 (850) 308 6742

WHATSAPP : +90 (533) 636 6742

Address : İncilipınar Mah. 36006 Cad.Kepkep Plaza
 No: 20/29 Kat: 4 Şehitkamil Gaziantep / TURKEY